

GT-521 PARTICLE COUNTER

SERVICE MAINTENANCE, REPAIR, SPARE PARTS

Met One Instruments, Inc.; 1600 NW Washington Blvd.; Grants Pass, OR. 97526
Ph. 541-471-7111; Fax 541-471-7116; E-mail: service@metone.com

REPAIR AND FIELD MAINTENANCE

This instrument must be returned to the factory for Calibration. The turn around time is generally less than five (5) working days.

RETURN AUTHORIZATION

Prior to any product return, please contact the factory for a return authorization number.

BASIC SERVICE(S)

ORDER #	DESCRIPTION	LIST PRICE
SVC0103	CALIBRATION, GT-521, GT-521S HAND HELD Includes Calibration and NIST Traceable Certificate Inspect Filter, Inspect, Clean Optical Sensor, Pump and Tubing. Cycle and Test Battery; Install any applicable upgrades. <i>Includes "As Received" and "As Left" NIST Cal certificates</i>	\$ 275.00
SVC3002	BENCH CAL, G3120 TEMP/RH SENSOR, NIST	\$ 65.00
SVC0109	ISO 21501-4 CALIBRATION	\$ 650.00

Parts Marked with an Asterisk (*) ahead of the part number will be replaced as needed during the calibration

REPLACEABLE PARTS

PART #	DESCRIPTION	LIST PRICE
3228	RS232 SERIAL CABLE	\$ 40.00
400100	CABLE, AC POWER	\$ 15.00
510640	POWER SUPPLY	\$ 75.00
G3110	ISO-KINETIC SAMPLING PROBE	\$ 55.00
8517	CARRYING CASE	\$ 45.00
993020	MINIATURE SCREW DRIVER	\$ 4.00
3305	GT-COMM SOFTWARE	\$ 100.00
GT-COM-9800	GT-COMM OPERATION MANUAL	\$ 25.00
GT-521-9800	GT-521 OPERATION MANUAL	\$ 25.00
G3111	ZERO FILTER	\$ 45.00
G3120	TEMPERATURE / HUMIDITY SENSOR	\$ 350.00
G3115	PRINTER FOR USE WITH HANDHELD MONITORS	\$ 450.00
750509	PRINTER PAPER	\$ 3.50
* 8538	IN-LINE FILTER ASSEMBLY	\$ 30.00
* 9285	Battery Pack, GT-521, Aerocet 531, GT-526	\$ 90.00
* 8650	PUMP ASSY	\$ 250.00