

BT-620 BENCH TOP PARTICLE MONITOR 1 CFM

SERVICE MAINTENANCE, REPAIR, SPARE PARTS

Met One Instruments, Inc.; 1600 NW Washington Blvd.; Grants Pass, OR. 97526
Ph. 541-471-7111; Fax 541-471-7116; E-mail: service@metone.com

REPAIR AND FIELD MAINTENANCE

This instrument must be returned to the factory for Calibration. The turn around time is generally less than five (5) working days.

RETURN AUTHORIZATION

Prior to any product return, please contact the factory for a return authorization number.

BASIC SERVICE(S)

ORDER #	DESCRIPTION
SVC0110	CALIBRATION, BT-620 BENCH TOP PARTICLE COUNTER Includes Calibration and NIST Traceable Certificate Inspect Filter, Inspect, Clean Optical Sensor, Pump and Tubing. Cycle and Test Battery; Install any applicable upgrades. <i>Includes "As Received" and "As Left" NIST Cal certificates</i>
SVC0109	ISO 21501-4 CALIBRATION

Parts Marked with an Asterisk (*) ahead of the part number will be replaced as needed during the calibration

REPLACEABLE PARTS

PART #	DESCRIPTION
81770	BATTERY CHARGER CABLE, BT-620
500784	USB CABLE
550065	RS232 CABLE, VARIOUS PRODUCTS.
81769	REPLACEMENT BATTERY PACK, BT-620
81751	BATTERY CHARGER, BT-620
81752	ISOKINETIC PROBE, BT-620
81754	ZERO FILTER, BT-620
G3120	TEMPERATURE / HUMIDITY SENSOR
750514	THERMAL PRINTER PAPER
81755	FLOW METER, BT-620