

# 804 PARTICLE COUNTER

## SERVICE MAINTENANCE, REPAIR, SPARE PARTS

Met One Instruments, Inc.; 1600 NW Washington Blvd.; Grants Pass, OR. 97526  
Ph. 541-471-7111; Fax 541-471-7116; E-mail: service@metone.com

### REPAIR AND FIELD MAINTENANCE

This instrument must be returned to the factory for Calibration. The turn around time is generally less than five (5) working days.

### RETURN AUTHORIZATION

Prior to any product return, please contact the factory for a return authorization number.

### BASIC SERVICE(S)

ORDER #	DESCRIPTION
SVC0039	CALIBRATION, 804 PARTICLE COUNTER Includes "As Received" and "As Left" NIST Cal certificates Clean & Inspect flow system. Check battery & charge system.
SVC0109	ISO 21501-4 CALIBRATION

Parts Marked with an Asterisk (\*) ahead of the part number will be replaced as needed during the calibration

### REPLACEABLE PARTS

PART #	DESCRIPTION
80530	FLOW METER KIT
390031	BATTERY CHARGER, 8.4 VDC, 1.5A, LITHIUM ION, SMART
400113	POWER CORD, GT-526S, AEROCET 531S
500787	USB CABLE, MALE A TO MINI MALE B
80450	PROTECTIVE BOOT
8517	CARRYING CASE
80419	INLET SCREEN KIT, 804
80846	ZERO FILTER KIT

