

GT-526S PARTICLE COUNTER

SERVICE MAINTENANCE, REPAIR, SPARE PARTS

Met One Instruments, Inc.; 1600 NW Washington Blvd.; Grants Pass, OR. 97526
Ph. 541-471-7111; Fax 541-471-7116; E-mail: service@metone.com

REPAIR AND FIELD MAINTENANCE

This instrument must be returned to the factory for Calibration. The turn around time is generally less than five (5) working days.

RETURN AUTHORIZATION

Prior to any product return, please contact the factory for a return authorization number.

BASIC SERVICE(S)

ORDER #	DESCRIPTION
SVC0106	CALIBRATION, GT-526 PARTICLE COUNTER Includes Calibration and NIST Traceable Certificate Inspect Filter, Inspect, Clean Optical Sensor, Pump and Tubing. Cycle and Test Battery; Install any applicable upgrades. <i>Includes "As Received" data and "As Left" data certificates</i>
SVC3002	BENCH CAL, G3120 TEMP/RH SENSOR, NIST
SVC0109	ISO 21501-4 CALIBRATION

REPLACEABLE PARTS

PART #	DESCRIPTION
3228	RS232 SERIAL CABLE
400113	POWER CORD, GT-526S, AEROCET 531S
500787	USB CABLE, MALE A TO MINI MALE B
G3110	ISO-KINETIC SAMPLING PROBE
8517	CARRYING CASE
390031	BATTERY CHARGER, 8.4 VDC, 1.5A, LITHIUM ION, SMART
80248	COMET TERMINAL SOFTWARE CD
9801	FLOW METER
G3111	ZERO FILTER
G3120	TEMPERATURE / HUMIDITY SENSOR
G3115	PRINTER FOR USE WITH HANDHELD MONITORS
750509	PRINTER PAPER
* 80458-1	BATTERY PACK, 7.4 VDC 2000 mAh
* 8650-2	PUMP ASSY., FURGUT
580336	INLINE FILTER, (Between Sensor and Pump)
* 8538	IN-LINE FILTER ASSEMBLY

GT-526S Standard Equipment



GT-526S
P/N GT-526S



Manual
P/N GT-526S-9800



Comet Software CD
P/N 80248



Calibration Certificate
P/N GT-526S-9600



Battery Charger
P/N 390031



Power Cord
P/N 400113



ISO-Kinetic Probe
P/N G3110



Zero Count Filter
P/N G3111



Carrying Case
P/N 8517



USB Cable
P/N 500787



Custom Serial Cable
P/N 3228


Printer



MOI P/N: G3115

GT-526 Optional Equipment

Temperature & Humidity Probe



MOI P/N: G3120

Dwyer™ Ball Flow Meter



MOI P/N: 9801