

# THE DOCK

# **DOCKING STATION**



The Dock, Met One Instruments, Inc.'s ALL-NEW docking station is a power and communication accessory for select handheld devices. Compatible with the new generation of MOI handheld particle counters, including the Aerocet 532 and DR-528, The Dock will charge the handheld unit's battery and provide RS-485, Ethernet, and WiFi communication options simultaneously.

Data transfer is faster and easier than ever with this sleek, all-aluminum unit, sure to be a feature on any desk or tabletop. The Dock rapidly charges new generation handheld devices without blocking the large feature display, allowing the user to view readings on-screen from a distance while the instrument is charging. A removable spacer enables the user to dock their handheld without having to remove the rubber boot.

### **SPECIFICATIONS**

## **THE DOCK HANDHELD DEVICE DOCKING STATION**

#### Power

Input Power:

8.4 V, 1.3 A

#### Physical

Weight: Size: 1 lb, 5.4 oz (0.6 kg) 3.125 in (7.9 cm) height x 4.525 in (11.5 cm) width x 3.5 in (8.9 cm) depth

# **COMPATIBLE UNITS:**

AEROCET 532 Handheld
Particle Monitor



DR-528 Handheld Particle
Counter





Specifications are subject to change at any time.